

Metabolic Labeling of rRNA with $^{32}\text{P}_i$

Modified from Pestov et al 2008, Current Protocols in cell biology

1. Plate cells in 6-well dish 1 day before experiment. For U2OS cells, plate 3×10^5 cells.
2. One hour before starting labeling reaction, replace cell growth media with 1.0 ml of prewarmed phosphate-free media (ThermoFisher Gibco 11971025) with dialyzed FBS (ThermoFisher Gibco 26400044).
3. Typically, if drug treating, add drugs to phosphate free media after 50 minutes of incubation.
4. Add $10 \mu\text{Ci}/\text{ml}$ of $^{32}\text{P}_i$ to media and incubate for 1 hour. Source radiation from Perkin-Elmer (NEX053H002MC) is $10 \mu\text{Ci}/\mu\text{l}$, thus, to the 1 mL of media, add $1 \mu\text{l}$ of $^{32}\text{P}_i$.
5. After 1 hour, remove $^{32}\text{P}_i$ containing media.
6. Wash 1x with prewarmed DMEM.
7. Add prewarmed DMEM. Make sure to add drugs if you intend on monitoring drug treatment.
8. At desired time-points, remove media, wash 1x with PBS
9. Add 0.5 mL of Trizol (Invitrogen) and extract with $100 \mu\text{l}$ of chloroform
10. Ethanol precipitate
11. Resuspend in $33 \mu\text{l}$ of RNA loading dye.
12. Run $10 \mu\text{l}$ on 1.2 % HEPES-EDTA-Formaldehyde gel overnight at 65V.
13. Dry gel and expose to film.

10X HEPES-EDTA Buffer

200 mM HEPES
10 mM EDTA
Adjust pH to 7.8

RNA loading buffer (make fresh)

Deionized formamide (95%)	7 μl
Formaldehyde (40%)	2 μl
10X H-E Buffer	1 μl
EtBr (400 $\mu\text{g}/\text{ml}$)	1 μl
0.5% Bromophenol Blue	1 μl